

# The **FastPack<sup>®</sup>** IP System

## PRODUCT SPECIFICATIONS

TECHNOLOGY:  
Paramagnetic Particle  
Chemiluminescence Immunoassay

OPERATING TEMPERATURE:  
15°C (59°F) to 32°C (90°F)

OPERATING HUMIDITY:  
10% to 80% Relative Humidity

DISPLAY OPTION:  
2-line Alphanumeric Display

DATA OUTPUT:  
RS 232, 9600 Baud, ASCII Text

FIELD SERVICE OPTIONS:  
ROM Card Firmware Upgrade

ELECTRICAL REQUIREMENTS:  
100-240 ± 10% VAC, 1.9 A

POWER CONSUMPTION:  
100 Watts (Maximum)

SIZE:  
13"W x 9"H x 12"D (33 x 23 x 30 cm)

WEIGHT:  
27 lbs (12.25 kg)



**Qualigen<sup>®</sup>**